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Instructions for use

Matrix Flow 70 the transparent addition-curing silicone out of a cartridge. Ideally for the silicon matrix technique, duplication technique and the flask technique in the anaxform flasks. The high shore hardness of 70 is ideal for the transfer of a wax-up directly into the mouth. Its easy handling, good flowability, and the high transparency ensure a quick polymerization of light-curing materials.

Indication:

- silicon matrix and die for light-curing composite
- transfer of a Wax-up in Mock-up
- cendric registration

Application:

To guarantee an optimum mixing proportion, make sure before attaching the mixingcanula, that both components (base and catalyst) steadily stamp out and spend a tiny amount on a mixing block. Please repeat this before every use. Afterwards put the mixingcanula on the cartouche. Before using the matrix Flow 70 clear please make sure that the contact surfaces are clean, grease-free and dry. During the application make sure that the tip of the mixingcanula is dipped in the material to avoid the encircling of air bubbles

Technical data:

Mixing ratio	1 : 1
Shore hardness after one hour	72
Working time at 23 °C / 73,4 °F	ca. 1,5 min.
Curing time at 23 °C / 73,4 °F	ca. 6 min.
C. time under pressure 2,5 bar	ca. 6 min.
Volume mass at 23 °C / 73,4 °F	ca. 1,55
linear dimension expending	< 0,185 %
Deformation under pressure % ISO	0,1%

Storage:

Keep material closely, separate and store them in a ventilated place at 5-27°C (41-77°F).

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anax dent GmbH Olgastr. 120a D-70180 Stuttgart