

# Bearbeitungssets / Tool kits



**MODELLHERSTELLUNG + GIPS BEARBEITUNGSSET**  
**MODEL PREPARATION + PLASTER TOOL-KIT**

300187

Order-Nr.

### Inhalt/Content

- 01 80055 5'-10000 Upm/rpm
- 02 153520 5'-10000 Upm/rpm
- 03 60017 5'-15000 Upm/rpm
- 04 509101 15'-25000 Upm/rpm
- 05 60024 3'-15000 Upm/rpm
- 06 60028 5'-10000 Upm/rpm
- 07 60029 5'-10000 Upm/rpm



**TITAN BEARBEITUNGSSET**  
**TITANIUM TOOL-KIT**

300117

Order-Nr.

- 01 708001 15'-30000 Upm/rpm
- 02 700801 15'-30000 Upm/rpm
- 03 702602 15'-30000 Upm/rpm
- 04 704602 15'-30000 Upm/rpm
- 05 704802 15'-30000 Upm/rpm
- 06 702402 15'-30000 Upm/rpm
- 07 702102 15'-30000 Upm/rpm
- 08 702202 15'-30000 Upm/rpm
- 09 704902 15'-30000 Upm/rpm
- 10 80220 3'- 5000 Upm/rpm
- 11 80222 3'- 5000 Upm/rpm
- 12 80305 1'- 3000 Upm/rpm
- 13 80301 1'- 3000 Upm/rpm



**NEM (K+B) BEARBEITUNGSSET**  
**NON PRECIOUS METALS (CROWN-BRIDGE) TOOL-KIT**


300157

Order-Nr.

- 01 60052 10'-20000 Upm/rpm
- 02 133520 3'-15000 Upm/rpm
- 03 201801 15'-30000 Upm/rpm
- 04 201402 15'-30000 Upm/rpm
- 05 201602 15'-30000 Upm/rpm
- 06 201202 15'-30000 Upm/rpm
- 07 118050 15'-30000 Upm/rpm
- 08 304802 20'-30000 Upm/rpm
- 09 80014 5'-20000 Upm/rpm
- 10 117014 15'-30000 Upm/rpm
- 11 54217 15'-30000 Upm/rpm
- 12 304402 15'-30000 Upm/rpm
- 13 305002 15'-30000 Upm/rpm
- 14 60032 3'-12000 Upm/rpm
- 15 80305 1'-3000 Upm/rpm
- 16 80301 1'-3000 Upm/rpm

15001 Reinigungsstein / cleaning stone

# Bearbeitungssets / Tool kits



**KERAMIK BEARBEITUNGSSET  
CERAMICS TOOL-KIT**

300137



Order-Nr.

## Inhalt/Content

01	60042608	3'-30000	Upm/rpm
02	60110258	3'-30000	Upm/rpm
03	60199258	3'-30000	Upm/rpm
04	60250168	3'-30000	Upm/rpm
05	60166188	3'-30000	Upm/rpm
06	60225188	3'-30000	Upm/rpm
07	60001238	3'-30000	Upm/rpm
08	60001148	3'-30000	Upm/rpm
09	60014238	3'-30000	Upm/rpm
10	112221	5'-10000	Upm/rpm
11	133622	5'-10000	Upm/rpm
12	812050	5'-10000	Upm/rpm
13	800120	3'-12000	Upm/rpm
14	801911	2'-3000	Upm/rpm
15	801920	2'-3000	Upm/rpm
16	80041	2'-3000	Upm/rpm



**EDELMETALL (K+B) BEARBEITUNGSSET  
PRECIOUS METALS (CROWN+BRIDGE) TOOL-KIT**

300147



Order-Nr.

01	302801	15'-30000	Upm/rpm
02	305002	15'-30000	Upm/rpm
03	302602	15'-30000	Upm/rpm
04	304602	15'-30000	Upm/rpm
05	704802	15'-30000	Upm/rpm
06	117014	15'-30000	Upm/rpm
07	117006	15'-30000	Upm/rpm
08	118050	15'-30000	Upm/rpm
09	117211	15'-30000	Upm/rpm
10	272125	3'-5000	Upm/rpm
11	80222	3'-5000	Upm/rpm
12	80301	1'-3000	Upm/rpm
13	812050	5'-10000	Upm/rpm
14	80305	1'-3000	Upm/rpm



**KOMPOSIT + KUNSTSTOFFVERBLENDUNGSSET  
TOUGH ACRYLIC VENEERS + COMPOSITES TOOL-KIT**

300227



Order-Nr.

01	700801	15'-30000	Upm/rpm
02	704901	15'-30000	Upm/rpm
03	702602	15'-30000	Upm/rpm
04	162221	3'-5000	Upm/rpm
05	118080	15'-30000	Upm/rpm
06	60019	5'-10000	Upm/rpm
07	801911	2'-3000	Upm/rpm
08	801940	2'-3000	Upm/rpm
09	801960	2'-3000	Upm/rpm
10	803060	2'-3000	Upm/rpm
11	80301	1'-3000	Upm/rpm
12	80041	1'-3000	Upm/rpm
13	600501	5'-10000	Upm/rpm
15001	Reinigungsstein / cleaning stone		

# Bearbeitungssets / Tool kits



**KUNSTSTOFF BEARBEITUNGSSET  
ACRYLICS TOOL-KIT**

300127



Order-Nr.

## Inhalt/Content

01	509001	15'-30000 Upm/rpm
02	407001	15'-30000 Upm/rpm
03	406602	15'-30000 Upm/rpm
04	117006	15'-30000 Upm/rpm
05	117023	15'-30000 Upm/rpm
06	118050	15'-30000 Upm/rpm
07	80040	5'-10000 Upm/rpm
08	80028	5'-10000 Upm/rpm
09	80225	3'-5000 Upm/rpm
10	152211	3'-5000 Upm/rpm



**FRÄSTECHNIK SET 0°  
MILLING TECHNIQUE KIT 0°**

300247



Order-Nr.

01	115712L	3'-6000 Upm/rpm
02	115422L	6'-12000 Upm/rpm
03	115602L	6'-12000 Upm/rpm
04	115802L	6'-12000 Upm/rpm
05	80071	5'-10000 Upm/rpm
06	80072	5'-10000 Upm/rpm
07	80073	5'-10000 Upm/rpm
08	115722L	3'-6000 Upm/rpm
09	115432L	6'-12000 Upm/rpm
10	115612L	6'-12000 Upm/rpm
11	115812L	6'-12000 Upm/rpm
12	115732L	3'-6000 Upm/rpm
13	115442L	6'-12000 Upm/rpm
14	115622L	6'-12000 Upm/rpm
15	115822L	6'-12000 Upm/rpm
16	80061	5'-10000 Upm/rpm
17	80062	5'-10000 Upm/rpm
18	80063	5'-10000 Upm/rpm



**KONOMETRIESET 2° und 4°  
MILLING TECHNIQUE KIT 2° und 4°**

300248



Order-Nr.

01	115727L	3'-6000 Upm/rpm
02	115472L	6'-12000 Upm/rpm
03	115627L	6'-12000 Upm/rpm
04	115137L	6'-12000 Upm/rpm
05	80461	5'-10000 Upm/rpm
06	80462	5'-10000 Upm/rpm
07	80463	5'-10000 Upm/rpm
08	115717L	3'-6000 Upm/rpm
09	115462L	6'-12000 Upm/rpm
10	115617L	6'-12000 Upm/rpm
11	115127L	6'-12000 Upm/rpm
12	80261	5'-10000 Upm/rpm
13	80262	5'-10000 Upm/rpm
14	80263	5'-10000 Upm/rpm