


Gipstrennscheiben / Plaster discs

Anwendung/Application



- Gips
Plaster
- Kunstharze
Synthetic resins

Körnung / grain size
ISO 514
ISO-Nr.
806 104 ...
Figur-Nr. / Shape #
ø 1/10 mm
VE / packaging unit

fein, roter Ring / fine, red ring

MINI	PLASTER DISC	MAXI
368	5' - 10000 Upm/rpm (low)	368
320		400
1		1


153222

154022

Order-Nr.



- Gips
Plaster
- Kunstharze
Synthetic resins



fein, roter Ring / fine, red ring

HIGH-TECH	PLASTRO	HIGH-TECH
379	5' - 10000 Upm/rpm (low)	379
350		450
1		1

153522

153722


154522

Order-Nr.



Gipstrennscheiben / Plaster discs

Galvanische
Diamantierung
*Electroplated
diamonds*



Körnung / grain size
ISO 504
ISO-Nr.
806 000 ...
Figur-Nr. / Shape #
ø 1/10 mm
VE / packaging unit

superfein / superfine

GIPSTRENNSCHEIBE / PLASTER CUT
3' - 10000 Upm/rpm (high)

Order-Nr.

368	368	368	368	368
750	750	750	750	800
1	1	1	1	1
1575142	1575162	1575202	1575222	1680162
x = 14mm H7	x = 16mm H7	x = 20mm H7	x = 22mm H7	x = 16mm H7

Circular gesinterte Galvano Scheibe
Circular sintered long life Galvano disc



Körnung / grain size
ISO 514
ISO-Nr.
806 000 ...
Figur-Nr. / Shape #
ø 1/10 mm
VE / packaging unit

fein / fine

PLASTER-CUT LONGTIME
3' - 10000 Upm/rpm (high)

Order-Nr.

1580162